

## Tissue Microarray - Human Cervical Carcinoma

Cat.-No.: 401 2203 Sample Datasheet

Slide Label						
	a	b	c	d	e	f
1	●	●	●	●	●	●
2	●	●	●	●	●	●
3	●	●	●	●	●	●
4	●	●	●	●	●	●
5	●	●	●	●	●	●
6	●	●	●	●	●	●
7	●	●	●	●	●	●
8						
9	●	●	●	●	●	
10	●	●	●	●	●	

Technical Information: 52 spots

- Spot diameter: 2.0 mm
- Fixation in 4% neutral buffered formaldehyde solution
- Paraffin embedded

Tissue type validated by immunohistochemistry

**In vitro laboratory use only.**

Not intended for any human or animal diagnostic or therapeutic use.

Position	Diagnostics	Sex	Age	Stage pT	Stage pN	Stage pM	Grade
1a	SCC	f	48	4	0	X	3
1b	SCC	f	30	2b	0	X	3
1c	SCC	f	31	X	1	1	3
1d	SCC	f	39	1b	0	X	3
1e	ADC	f	46	2b	0	X	2
1f	SCC	f	32	4	0	X	3
2a	SCC	f	69	2	X	X	2
2b	SCC	f	45	1b	0	X	2
2c	SCC	f	61	1b	0	X	2
2d	SCC	f	41	1b	0	X	2
2e	SCC	f	44	1b	0	X	2
2f	SCC	f	52	2b	1	X	3
3a	SCC	f	43	1b	0	X	2
3b	SCC	f	43	2b	0	X	2
3c	SCC	f	38	y1b	0	X	3
3d	SCC	f	36	2b	1	X	2
3e	SCC	f	62	2	0	X	2
3f	SCC	f	25	1b	1	X	3
4a	SCC	f	47	2b	1	1	2
4b	lymph node metastasis	f	47	2b	1	1	2
4c	SCC	f	72	2b	X	X	3
4d	SCC	f	40	1b	0	X	2
4e	SCC	f	38	1 a1	X	x	2
4f	SCC	f	68	1 b1	0	x	2
5a	SCC	f	50	2b	1	X	3
5b	lymph node metastasis	f	50	2b	1	X	3
5c	SCC	f	43	1b	1	X	2

**In vitro laboratory use only.**

Not intended for any human or animal diagnostic or therapeutic use.

Position	Diagnostics	Sex	Age	Stage pT	Stage pN	Stage pM	Grade
5d	SCC	f	69	1	X	X	3
5e	SCC	f	56	2b	0	X	3
5f	SCC	f	32	2b	1	X	2
6a	SCC	f	46	2a	0	X	3
6b	SCC	f	70	2b	0	X	2
6c	SCC	f	62	X	X	1	3
6d	Peritumorous inflammatory tissue	f	49	y2b	X	1	3
6e	SCC	f	35	1b2	0	X	2
6f	SCC	f	69	4	0	X	2
7a	SCC	f	42	2b	1	X	2
7b	SCC	f	58	2b	X	X	3
7c	SCC	f	44	1b1	X	X	2
7d	SCC	f	33	1b1	0	X	2
7e	SCC	f	37	1b1	0	X	3
7f	SCC	f	59	1bi	0	X	3
9a	CIN III lesion	f	50				
9b	CIN III lesion	f	63				
9c	Cervix mucosa near CIN III lesion	f	40				
9d	CIN III lesion	f	29				
9e	Cervix mucosa near CIN III lesion	f	26				
10a	normal mucosa	f	30				
10b	normal mucosa	f	32				
10c	normal mucosa	f	35				
10d	normal cervix, fibrous tissue	f	36				
10e	normal mucosa	f	35				

**In vitro laboratory use only.**

Not intended for any human or animal diagnostic or therapeutic use.