

## Tissue Microarray - Human Autoimmune Tissue

Cat.-No.: 401 3301 Sample Datasheet

Slide Label						
	a	b	c	d	e	f
1	●	●	●	●	●	●
2	●	●	●	●	●	●
3	●	●	●	●		
4	●	●	●	●		
5	●	●	●	●		
6	●	●	●			
7	●	●	●	●	●	●
8	●	●	●	●	●	●
9	●	●	●	●		

Technical Information: 43 spots

- Spot diameter: 2.0 mm
- Fixation in 4% neutral buffered formaldehyde solution
- Paraffin embedded

Tissue type validated by immunohistochemistry


**In vitro laboratory use only.**

Not intended for any human or animal diagnostic or therapeutic use.

Position	Diagnosis	Tissue	Sex	Age
1a	Hashimoto-Thyreoiditis	Thyroid	f	37
1b	Hashimoto-Thyreoiditis	Thyroid	f	37
1c	Hashimoto-Thyreoiditis	Thyroid	f	13
1d	Hashimoto-Thyreoiditis	Thyroid	f	13
1e	Hashimoto-Thyreoiditis	Thyroid	f	49
1f	Hashimoto-Thyreoiditis	Thyroid	f	49
2a	Primäre biliäre Zirrhose	Liver	m	47
2b	Primäre biliäre Zirrhose	Liver	m	47
2c	Sjögren	Submandibularis	f	72
2d	Sjögren	Submandibularis	f	72
2e	Sjögren	Labium oris	f	47
2f	Sjögren	Labium oris	f	47
3a	Sinusitis & Eosinophilie	Nasal mucosa	m	62
3b	Sinusitis & Eosinophilie	Nasal mucosa	m	62
3c	Sinusitis & Eosinophilie	Nasal mucosa	m	62
3d	Sinusitis & Eosinophilie	Nasal mucosa	m	62
4a	Rheumatoide Arthritis	Synovial	f	59
4b	Rheumatoide Arthritis	Synovial	f	59
4c	Psoriasis	Synovial	m	46
4d	Psoriasis	Synovial	m	46
5a	Morbus Crohn	Sigma (muscularis)	f	17
5b	Morbus Crohn	Sigma (mucosa)	f	17
5c	Sarkoidose	Lung	f	67
5d	Sarkoidose	Lung	f	67
6a	Granulomat. Thyreoiditis de Quervain	Thyroid	f	52
6b	Granulomat. Thyreoiditis de Quervain	Thyroid	f	52

**In vitro laboratory use only.**

Not intended for any human or animal diagnostic or therapeutic use.

 Sample Datasheet  
Cat.-No.: 401 3301

Position	Diagnosis	Tissue	Sex	Age
6c	Lupus erythematodes	Dermis	f	65
7a	normal	Thyroid	m	53
7b	normal	Thyroid	m	53
7c	normal	Liver	m	61
7d	normal	Liver	m	61
7e	normal	Salivary gland	m	40
7f	normal	Salivary gland	m	40
8a	normal	Dermis	m	41
8b	normal	Dermis	m	41
8c	normal ( Sinusitis w/o. Eosinophilie)	Nasal mucosa	m	36
8d	normal ( Sinusitis w/o. Eosinophilie)	Nasal mucosa	m	36
8e	normal	Sigma	m	58
8f	normal	Sigma	m	58
9a	normal	Lung	m	47
9b	normal	Lung	m	47
9c	normal	Synovial	m	30
9d	normal	Synovial	m	30

**In vitro laboratory use only.**

Not intended for any human or animal diagnostic or therapeutic use.