

Tissue Microarray - Human Synovitis Tissue

Cat.-No.: 401 3201 Sample Datasheet

Slide Label						
	a	b	c	d	e	f
1	●	●	●	●	●	●
2	●	●	●	●		
3	●	●	●	●	●	●
4	●	●	●	●	●	●
5	●	●	●	●	●	●

Technical Information: 28 spots

- Spot diameter: 2.0 mm
- Fixation in 4% neutral buffered formaldehyde solution
- Paraffin embedded

Tissue type validated by immunohistochemistry /classified according to synovitis score by Krenn (Pathol Res Pract. 2002)

In vitro laboratory use only.

Not intended for any human or animal diagnostic or therapeutic use.

📖 Sample Datasheet
Cat.-No.: 401 3201

Position	Tissue	Diagnosis	Score	Sex	Age
1a	wrist	RA	6 / 9	f	65
1b	wrist	RA	6 / 9	f	65
1c	ankle	RA	6 / 9	m	70
1d	ankle	RA	6 / 9	m	70
1e	shoulder	RA	8 / 9	f	32
1f	shoulder	RA	8 / 9	f	32
2a	knee	RA	8 / 9	m	62
2b	knee	RA	8 / 9	m	62
2c	joint	RA	7 / 9	f	59
2d	joint	RA	7 / 9	f	59
3a	knee	PSA	5 / 9	m	45
3b	knee	PSA	5 / 9	m	45
3c	knee	normal tissue	0 / 9	m	56
3d	knee	normal tissue	0 / 9	m	56
3e	SCC III	normal tissue	1 / 9	f	79
3f	SCC III	normal tissue	1 / 9	f	79
4a	hip	OA	5 / 9	f	76
4b	hip	OA	5 / 9	f	76
4c	joint	OA	4 / 9	f	88
4d	joint	OA	4 / 9	f	88
4e	knee	OA	3 / 9	m	64
4f	knee	OA	3 / 9	m	64
5a	hip	OA	6 / 9	f	77
5b	hip	OA	6 / 9	f	77
5c	hip	OA	4 / 9	f	73
5d	hip	OA	4 / 9	f	73
5e	joint	OA	3 / 9	f	63
5f	joint	OA	3 / 9	f	63

In vitro laboratory use only.

Not intended for any human or animal diagnostic or therapeutic use.